



TiVi700 Analyzer Main Manual

User Manual 1.5 February 2019

TISSUE VIABILITY IMAGER

Information in this document is subject to change without notice and does not represent a commitment on the part of WheelsBridge AB. No part of this manual may be reproduced or transmitted in any form or by any means, electronically or mechanically, including photocopying and recording, for any purpose without the express written permission of WheelsBridge AB.

Any use or application of WheelsBridge AB products must be performed by properly trained and qualified personnel. WheelsBridge AB takes no responsibility for action taken on the basis of information derived from the use of the product described in this manual or from consequential damages in connection with the furnishing, performance or use of this material.

Licensee's limited rights are governed by the license agreement between WheelsBridge AB and MathWorks. Licensee may not modify or remove any license agreement file (MarthWorks or third party) that is included with the MCR Libraries ("MCR Library License"). Licensee of this Application accepts the terms of the MCR Library License.

The *TiVi700 Tissue Viability Imager* is patented and/or patents are pending in several countries.

The *TiVi700 Tissue Viability Imager* is not registered as a Medical Device. It is intended for research applications only.

WheelsBridge is a registered trade mark.

CONTENT

Part I	Introduction
Part 2	Menu Items
Part 3	Photo frame
Part 4	Image frame
Part 5	Select Photos Frame
Part 6	Step Through Frame
Part 7	Crop Images frame
Part 8	Movie frame
Part 9	View Frame
Part 10	Single Image Window
Part 11	Group of Images Window
Part 12	Curve Tracker Window
Part 13	Gallery Window
Part 14	Image Slicer Window
Part 15	Cross Sectioner Window
Part 16	Profiler Window
Part 17	Wizard
Part 18	Chart Window
Part 19	Project Window
Part 20	Excel Files Generated
Part 21	Specifications
Part 22	Demo Assistant
Part 23	DropData
Part 24	Organizer

Part 25 Software Installation Procedure

Part 1 – Introduction

The **TiVi700 Analyzer** is the analysis part of the **TiVi700 Tissue Viability Imaging** system. The **TiVi700 Analyzer** analysis sequences of photos and creates a TiVi image mapping the skin red blood cells concentration for each photo. The **TiVi 700 Analyzer** is composed of two main pages – the **Show/Crop** page and the **Movie/View** page – described in further detail below.

The TiVi700 software runs on 64 Bit computers only.

Part 2 – Menu Items

- 1. The **File** menu.
 - a. Camera Opens the TiVi700 Camera software for capturing of TiVi photos. The TiVi700 Camera software is further described in the TiVi700 Camera manual.
 - b. **Open A Photo** opens a new **Photo** in the **Photo** window. This window can be used as a reference to the **Photo** and **Image** being displayed in the **TiVi700** main window. Multiple **Photo** windows can be kept open simultaneously, moved to outside the main window and rescaled to a suitable size. To rescale an open window place the mouse pointer at the border of the window (the pointer changes to a double arrow), press the left mouse button and drag the mouse to a new position.



c. **Open TiVi Chart** - opens a **TiVi Chart** saved from inside the TiVi **Chart** window.

- d. **Projects** → **New TiVi Project** opens a new **TiVi Project**.
- e. **Projects** → **Open TiVi Project** opens an existing **TiVi Project**.
- f. Save As for saving the actual sequence of Photos as displayed in the TiVi700 Analyzer window. This command is useful for saving the sequence of Photos in a different folder or under a new name. Saving regions of interest (ROIs) is done by use of the Crop & Save command in the Crop Images frame. All relevant parameter settings are saved together with the Photos in such a way that the next time the Photos are loaded, the associated TiVi Images will be displayed under the view angle and with the same mask settings, upper and lower limit as when they were stored. This is equivalent to saving the actual Mask settings of the Images.
- g. Load ROIs loads a selected set of ROIs. This set of ROIs is accepted only if the size of the actual **photo** is the same as the **photo** in which the **ROIs** were originally generated.
- h. Save ROIs saves the present set of ROIs to file.
- Photo resolution opens a new window that allows the user to save a sequence of the actual photos at lower resolution than the original photos. Working with the sequence of photos at lower resolution speeds up the data analysis procedures and is useful when a first overview of data is made. In the Photo Resolution window the user can set the desired resolution of the photos and save these photos to file under a new name.
- j. **Photo rotation** opens a new window that allows the user to rotate the photos a selected number of degrees and save a sequence of actual photos rotated to this degree. These rotated photos can be saved to file under a new name.
- k. **TiVi Pre-Processor** opens a new window that allows the user to align the object in the individual photo before analysis. For further information please see the separate **TiVi Pre-Processor** manual.
- Export Image exports the raw data associated with the Image displayed to a file that can be read by Excel® or by MATLAB® for further processing. Select if the entire sequence of photos or just the actual image are to be saved. Note that the Excel® decimal delimiter must be a dot and not a comma.
- m. Exit opens a dialog window that prompts the user to confirm the close of the TiVi700 Analyzer window.
- 2. The **Edit** menu.
 - a. Activate Pigm Comp allows the user to suppress the influence of pigmented spots on the TiVi-image. Fine-tuning is made by the **Pigmentation** Compensation on slider.
 - b. Deactivate Pigm Comp deactivates the pigmentation compensation settings.

- c. Copy Photo opens a Photo window from which the Photo and Photo Data can be copied to the Clipboard.
- d. Copy Image opens an Image window from which the Image and Image Data can be copied to the Clipboard.
- e. **Circle ROI / Rect ROI** sets the region of interest to circular or rectangular shape.
- f. Freehand ROI draws ROI by use of the mouse to user-selected shape.
- g. Common ROI all ROI operations are applied to all photos in the sequence.
- h. Individual ROI No Lock on when first drawn, the ROI is applied to all photos in the sequence. When a ROI is moved or stretched, this maneuvre applies only to the ROI in the actual Photo. By this arrangement a slight movement of the object in one or more photos can be compensated for by moving the ROI(s) to the correct position. When a ROIs is removed it is removed from all photos in the sequence. The individual ROIs are stored in a 3-dimensional matrix saving the actual photo file number, the ROI number and the ROI coordinates. When a ROI is saved, this matrix is saved with a prefix identical to the actual photo prefix to which a user selected sub-prefix is added.
- i. Individual ROI Lock on White Spot to keep the ROIs locked on specific positions on the target tissue even if the target is moving in relation to the camera. Start with placing a round white marker (a paper clip with a diameter of about 5 mm will do well) on the target tissue before capturing the photos. In the analysis with the setting to lock on to the white spot the system automatically finds the brightest spot and uses this as an anchor point for the ROIs in each photo. Draw ROIs and then click the Automatic button to scan through the photos. The ROIs should now appear locked in relation to the white spot in the photos. This feature is particularly useful when making ROI analysis on moving targets and assists in calculating the average ROI values from the same anatomical sites even though the target position may have changed in relation to the camera.
- j. Adjust ROIs sets the size of all ROIs to the size of the last set (white) ROI. The upper left corner position remains unchanged in all ROIs.
- k. Set ROI opens the ROI setting window from which it is possible to manually set the upper left and lower right coordinates of the last selected (white) ROI.
- **I. ROI Naming –** use to give every ROI a unique name.
- m. Remove Selected Area to remove areas selected by the ROIs from the photo. Draw ROIs around sections of the photo that is to be removed. Then click Remove Selected Area in the Edit pull-down menu. The new photos are saved in the folder RSA (Remove Selected Area) located in the folder of the

original photos. The new photos are given the prefix "RSA" followed by their original name. These new photos can be treated as standard TiVi-photos and when they are uploaded the selected areas are displayed in turquoise colour in the photo and as background in the corresponding TiVi image. This feature is useful for eliminating areas in the original photos that may distort the result such as a red spot characterizing an injection trauma at studies of the vaso-active response to an injected substance as wells as for elimination artificially coloured skin. This function operates also in cascade, i.e. the generated RSA photo can be used as the original photo to further remove new selected areas.

- 3. The **Page** Menu.
 - a. Crop Images to set the system to display the Step Through and the Crop Images frames.
 - b. Movie to set the system to display the Movie and the View frames.
- 4. The Analysis Menu
 - a. Single Image displays the Single Image window.
 - b. Group of Images displays the Group of Images window.
 - c. Curve Tracker displays the Curve Tracker window.
 - d. **Gallery** displays the **Gallery** window.
 - e. **Image Slicer** displays the **Image Slicer** window.
 - f. Cross Sectioner displays the Cross Sectioner window.
 - g. **Profiler** displays the **Profiler** window.
- 5. The Wizard Menu to display the Wizard for automatic analysis of TiVi data.
- 6. The DropData Menu
 - a. **Data to Sheet** to open a **DropData** window for inserting image data into a spread-sheet.
 - b. **Organizer** to open the Organizer table holding test subject data.
- 7. The **Maintenance** Menu.
 - a. Window Size changes the TiVi700 Analyzer window size between Small, Large and Wide Window (available only if the computer display is wide screen and high resolution).

- b. Generate Subfolders / No Subfolder Generation when many photos are generated as a result of image processing procedures, all photos can be stored in separate folders with the same name as the Prefix of the photo name. This feature helps to sort the individual photos and keep them separated. Click the Maintenance pull-down menu reading No Sub-folder generation to turn off this feature and click Generate Sub-folders to turn the feature on again.
- 8. The **Help** Menu.
 - a. **About** displays the **About** window with information about the software version and WheelsBridge AB.
 - b. License displays the software license agreement.
 - c. Manuals electronic copies of the Quick Reference Manual, the Main Manual and the Pre-Processor Manual.
 - i. Click e.g. the **Main Manual**. This manual will now open inside the **Abode PDF-Reader** (this assumes that the **Abode PDFReader** is installed). The **Abode PDFReader** can be down-loaded free of charge from <u>http://www.adobe.com</u>).
 - ii. Diminish the Adobe Acrobate Standard window so both this window and the TiVi700 Analyzer window are simultaneously visible on the screen.



- iii. If you for example want to retrieve information about how the **First Photo** button works, select Edit→ Find, print "First Image" in the edit box and click **Enter** on the computer keyboard.
- iv. To find the next "**First Photo**" position in the text, click the right arrow next to the edit box.

- v. Click "X" to hide the edit text box.
- d. TiVi700 Presentation an overview presentation of the TiVi technology.
- e. **Demo Assistant** a full explanation of the *Demo Assistant* features.
- d. **Product Sheet** displays detailed TiVi product sheet in pdf-format.
- e. WheelsBridge WEB a direct link to the WheelsBridge WEB-site.
- f. **ToolTips** –displays (ON) or hides (OFF) the tooltips associated with the different controls and displays.
- g. Patents lists TiVi patent numbers
- 9. The **Language** Menu to set the actual language. This information is stored between sessions and is automatically retrieved when a new session is initiated.
- 10. The **Tool Boxes** Menu to start individual toolboxes (optional).
- 11. The **Demo Assistant** Menu to start the interactive **Demo Assistance** guiding the User on how to use the different features of the **TiVi700 Analyzer**,

Part 3 - Photo frame

- 1. The **Photo** loaded from file or being the result of a **Show ROI** click button.
- 2. The Photo Size (in pixels) is displayed next to the Photo.
- 3. The Actual Photo Indicator indicates in blue color the relative number of the Photo in the sequence that is currently displayed.
- 4. **Date** indicates the date when the photo was captured.
- 5. Time First indicates the time when the first Photo of the sequence was captured.
- 6. **Time** indicates the time when the currently displayed **Photo** was captured as well as the **Elapsed Time** from when the **First Photo** was captured.
- 7. Time Last indicates the time when the last Photo of the sequence was captured.
- 8. **ROIs** displays the number of active regions of interest.
- 9. Information about **Polarization Mode**, **Distance** to object and **Zoom** mode as stored with the photo.

10. When a region of interest (ROI) is drawn the **Distance** (in pixels) between two diametric corners of the ROI is displayed below the **TiVi-value**.

Part 4 - Image frame

- 1. The **Image** calculated by the *TiVi700* algorithm or being the result of a **Show ROI** button click.
- 2. The Data Values for the displayed Image are listed to the right of the Image. Data Values representing the entire Image or a ROI are listed in blue color, while Data Values representing areas selected with the Upper Limit and Lower Limit sliders are listed in red color. Data Values representing a Subtraction Image are displayed in turquoise color, while Data Values representing a Subtraction Image with values within the Lower and Upper Limit are presented in gray color. Double presentation (x/y) for Mean indicates mean value including background (x) and only background (y). Double presentation (x/y) for Max and Min indicates the highest/lowest value (x) respectively the median value of the 10 highest/lowest values. The latter is used by the system in Mask calculations in order to avoid the influence of single extreme values.
- 3. Following selection made by the **Lower** and **Upper Limit** sliders, the limits (in percentage) are displayed above the **Image**.
- 4. The **TiVi-value** of the point at which the pointer is pointing in the **Photo** or **Image**, is continuously updated and displayed.
- 5. The median RGB value triplet for the actual photo displayed in the **Photo Frame**. For a sequence of photos, the corresponding RGB value triplets can be exported to the **Group Of Images** window where they can be displayed in a separate colour space window and exported to an Excel® document.
- 6. The **Color Bar** is displayed to the right of the **Image** and reflects the setting of the **Color Span**.

Part 5 - Select Photos Frame

- First Photo button to select the First Photo file of the sequence of Photos to be investigated from the files listed in the Open Dialog Box. The file name is displayed in the associated First Photo edit box. By default the Last Photo is automatically selected as the last photo in the consecutive sequence. The Last Photo of the sequence is displayed.
- 2. Last Photo button to select the Last Photo file of the sequence of Photos to be investigated from the files listed in the **Open Dialog Box**. The file name is displayed in the associated **Actual Photo** edit box.

- 3. **Step** to select the step in the **Photo** sequence to be loaded. If for example the **Photos** are stored in files with names F-1000, F-1001....F-1020, setting **Step** to 5 loads **Photo** data for the files F-1000, F-1005, F-1010, F-1015 and F-1020. The value is selected by repeatedly clicking the **Step** button.
- 4. **Reverse –** reverses the order of updating **Step** when the **Step** button is successively clicked. The button name toggles between **Reverse** and **Forward**.
- 5. Clear All clears the Photo and Image frames and prepares the system for accepting new Photos.
- 6. Show First Photo and Show Last Photo by clicking one of these radio-buttons, the first or last Photo and Image pair to be displayed can be selected.
- 7. Actual Photo box shows the name of the Photo file currently displayed. To directly display a specific Photo in a selected sequence, print its name in the Actual Photo edit box and press the Enter button on the computer keyboard.
- 8. When the system is in **Subtraction** mode, the name of the file representing the **Image** to be subtracted is displayed in the **Subtract Photo** edit box.

Part 6 - Step Through Frame

- 1. **Manual** by successively clicking this button, the **Photos** and **Images** are displayed in their respective frame, while the actual file name is shown in the **Actual Photo** edit box.
- 2. **Reverse -** changes the order in which the **Photos** and **Images** are displayed when successively clicking the **Manual** button. The button name toggles between **Reverse** and **Forward**.
- 3. Reset displays the First Photo and First Image.
- 4. **Delay (sec)** dragging the slider to a new position sets a new **Speed**, as displayed in the **Edit** box, at which the **Image** sequence following clicking **Automatic** is presented. Keeping the pointer over the leftmost or rightmost arrow in the slider, with the left mouse button down, continuously moves the slider while the **Edit** box is updated to a new **Speed** setting. Alternatively print a new value in the **Edit** box and press **Enter**.
- 5. Automatic clicking this button makes the system step through the **Photo** sequence at a speed set by **Delay.** Activating **Automatic** enables the **Pause** and **Abort** buttons.
- 6. **Pause** clicking this button while in **Automatic** mode, halts the **Step Through** procedure and waits for the button, now labeled **Continue**, to be clicked again.
- 7. Abort clicking this button while in Automatic mode, aborts the Step Through procedure.

- 8. **Comments** opens the word pad associated with the actual photo / image displayed for further editing. When the edit process is finished click the **Save & Close** button in the **Comments 700** window to save the edited text as meta data linked to the actual photo.
- 9. **Color Span** for setting the color scale of the displayed image, without changing the actual data values. The color scale end-points are indicated in blue (lowest value) and red (highest value).

Part 7 - Crop Images frame

- 1. Show ROI after having selected a Region of Interest (ROI), using the following procedure
 - Place the pointer within the **Photo** or the **Image** area in the upper left corner of the **ROI** to be generated (if the **ROI** is to be generated in the **Image** and this is not displayed in **Top View** mode (azimuth = 0 and elevation = 90) the User is prompted to place the system in **Top View** mode click **OK** in the **Warning** window to proceed),
 - Move the pointer to the lower right corner of the **ROI**, while keeping the left mouse button pressed,
 - Outline the **ROI** by releasing the left mouse button,

defines and outlines a ROI.

The position of any **ROI** can be changed by the following procedure:

When the pointer is moved over the upper left or lower right corner of a **ROI** the pointer changes shape from arrow to double arrow. When the pointer is indicated by a double arrow, the **ROI** shape can be changed by moving the pointer by use of the mouse with the left mouse button pressed. In **Common ROI** mode these changes are applied to all photos in the sequence. In **Individual ROI** mode, these changes are applied only to the actual photo.

When the pointer is moved over the centre of a **ROI**, the pointer changes shape from arrow to four-headed arrow. When the pointer is indicated by a four-headed arrow, the **ROI** can be moved to a new position (without changing its shape) by moving the pointer by use of the mouse with the left mouse button pressed. In **Common ROI** mode these changes are applied to all photos in the sequence. In **Individual ROI** mode, these changes are applied only to the actual photo.

When the pointer is positioned over the centre of a **ROI** as indicated by a fourheaded arrow, clicking the left mouse button causes the **Erase ROI** window to open. If the **Yes** button is clicked the actual **ROI** is erased in all photos and the remaining **ROIs** are re-numbered and re-drawn.

The actual number of **ROIs** is indicated by the **ROIs** indicator in the upper left corner of the **Photo Frame**. This number should always be the same as the number of **ROIs** displayed in the **Photo**. If not, click the **Refresh** button.

Note that a new **ROI** can also be outlined in the **Image**, but its border is always indicated in the **Photo**.

Switching between rectangular and circular **ROIs** is made by clicking **Circle ROI** or **Rect ROI** in the **Edit** pull-down menu.

The last drawn (white) **ROI** can now be displayed by clicking the **Show ROI** button. Again clicking the same button, now labeled **Undo Show ROI**, returns the system to its original mode. If new **ROIs** are successively generated the latest generated **ROI** is always displayed in white color while other **ROIs** are displayed in different colors so the **ROIs** can be easily identified. Clicking the **Refresh** button clears the **Photo** from **ROI** layouts. The **ROIs** can be drawn either in the **Photo** or in the **Image**.

By use of **Adjust ROI** in the **Edit** pull-down menu, the size of all ROIs can be made identical to the last **ROI** generated (white outline). The upper left corner of all **ROIs** is kept at the original position.

The coordinates of the upper left and lower right corner positions of the last **ROI** generated (white) can be set manually from the **ROI700setting** window opened by the **Set ROI** command in the **Edit** pull-down menu.

Provided at least one **ROI** exists, more **ROIs** of the same size can be automatically generated by placing the pointer at the upper left corner of the **ROI** to be generated and then clicking the left mouse button.

- 2. Crop & Save with the system in the Show ROI mode, clicking the Crop & Save button, opens the Save As dialog box. Select a folder (or alternatively generate a new folder) and a *file name* and click the Save button to save the sequence of (cropped) Photos. The saved sequence of cropped Photos adds the *file name* as a prefix to the old file name. All relevant parameter settings are saved together with the cropped Photos in such a way that the next time the (cropped) Photos are loaded, they will be displayed under the view angle they were stored. Furthermore Upper and Lower Limits are saved with the Photos so the actual threshold settings are retrieved. This is equivalent to saving the actual Mask settings of the Images.
- 3. Overlay after having selected upper and lower threshold by way of the Upper Limit and the Lower Limit slider, clicking the Overlay button, displays the actual reduced image areas superimposed on the photo, the latter now displayed in gray scale. When the Automatic button is clicked the selected image area is displayed superimposed throughout the entire photo sequence.

- 4. Subtract click this button to subtract the Actual Image from all the other Images in the sequence pixel by pixel. Clicking Subtract evokes the Subtracted Photo edit box displaying the name of the file representing the Image to be subtracted. To leave the Subtraction Mode again click the same button now named Undo Subtract.
- Upper Limit dragging this slider to a new position as displayed in the associated Edit box, sets the upper limit of Data Values to be included in the Image and displays this Image. Alternatively print a new value in the Edit box and press Enter.
- Lower Limit dragging this slider to a new position as displayed in the associated Edit box, sets the lower limit of Data Values to be included in the Image and displays this Image. Alternatively print a new value in the Edit box and press Enter.
- 7. Mask FiFr (Fixed Fraction) selects the fraction of the actual Image set by Lower Limit and Upper Limit and applies this Mask to all Images in the sequence. Fr (Fraction) selects the fraction set by Lower Limit and Upper Limit and applies this selection when setting the individual Masks to the Images. Fi% (Fixed Percentage) uses a Mask based on the values set by Lower Limit and Upper Limit in the actual Image and applies this Mask to all Images in the sequence. % uses individual Masks to all Image based on the Lower Limit and Upper Limit values. Level uses the set Lower and Upper Limits as absolute levels.
- 8. **Refresh** clears the **Photos** from **ROI** selections.
- 9. **Reset Mask** resets the actual mask.

Addendum: The **Crop Images** frame is accessible from the **Show/Crop Page**. The **ROI** selection, **Upper Limit** and **Lower Limit** selection may be combined to select **Image** data within a certain range and physical region of interest in the **Photo**. On the basis of this selection, the **Data Values** are calculated for further presentation in the **Group of Images** window. The **Mask** is particularly useful for automatically setting the **Region of Interest** and forms the core of increasing productivity in many applications. The possible settings (FiFr, Fr, Fi%,% and Level) is further illustrated in the Figure below. If an image contains single high values the Fi% and % setting may display a low number of active points within the **Mask** because very few values reach the set percentage **Limit** of the highest value. This problem is overcome when the FiFr and Fr settings are used because these settings will always display a number of points within the **Mask** corresponding to the fraction of the entire number of points in the image as set by the **Limit**.

Percentage(%)

Fraction (Fr)

	Fr	%
ixed	Mask calculated individually for each separate image	Mask calculated individually for each separate image
Non-F	Limit set as the fraction of all values	Limit set as percentage of highest value in the image
	FiFr	Fi%
ed	Mask always in the same place for all images as set in the reference image.	Mask always in the same place for all images as set in the reference image.
Fix	Limit set as the fraction of all values	Limit set as percentage of highest value in the image

Part 8 - Movie frame

- 1. Create Movie creates a Movie from the sequence of Images selected. While building the Movie, the system displays the different Image frames that will combine into the resulting Video Clip.
- 2. Interpol Video generates a video from the actual image sequence with a length (in seconds) that is set in the Video-Length window. Insert the desired length of the video in the edit text box and click the OK-button in the Video-Length window. The video is generated by calculating interpolated images in between the real images in the sequence. This procedure is particularly useful when a video is to be presented in e.g. a Power Point presentation.
- 3. Save Movie saves the created Movie in a file and folder set in the Save As dialog window.
- 4. Load Movie loads a previously created and saved Movie.
- 5. Speed dragging the slider to a new position sets a new Speed, as displayed in the Edit box, at which the Movie runs following clicking the Start Movie button. Keeping the pointer over the leftmost or rightmost arrow in the slider, with the left mouse button down, continuously moves the slider while the Edit box is updated to a new Speed setting. Alternatively print a new value in the Edit box and press Enter.
- 6. Start Movie –starts the Movie created by or loaded by the system.
- 7. Save As Images splits up a video-clip into individual photos and saves these to file.

Part 9 - View Frame

- 1. Azimuth dragging the slider to a new position changes the azimuth presentation angle as displayed in the associated Edit box. Keeping the pointer over the leftmost or rightmost arrow in the slider, with the left mouse button down, rotates the Image while the Edit box is continuously updated. The Image is rotated stepwise by successively clicking any of the arrows. Alternatively, a new value can be written in the Edit box. Activate the new settings by pressing Enter.
- 2. Elevation dragging the slider to a new position changes the elevation presentation angle as displayed in the associated Edit box. Keeping the pointer over the leftmost or rightmost arrow in the slider, with the left mouse button down, rotates the Image while the Edit box is continuously updated. The Image is tilted stepwise by successively clicking any of the arrows. Alternatively, a new value can be written in the Edit box. Activate the new settings by pressing Enter.
- 3. **Plot Density** dragging the slider to a new position changes the density of the plot grid associated with the **Image** presentation as displayed in the **Edit** box. Selecting a low **Plot Density**, speeds up the update procedure of the **Image** while at the same time the resolution is reduced, without affecting the **Data** values.
- 4. **Rotate Azim** clicking this button rotates the **Image** successively from 0 to 360 degrees. A low **Plot Density** speeds up the procedure. Interrupt the rotation procedure by clicking the button now labeled **Stop Rotate**.
- 5. **Rotate Elev** clicking this button changes the view angle under which the **Image** is seen over the range 0 to 90 degrees. A low **Plot Density** speeds up the procedure. Interrupt the rotation procedure by clicking the button now labeled **Stop Rotate**.
- 6. **Plot Method** selecting **Surface**, **Interpolated** and **Mesh** respectively displays the **Image** using different plot methods. The default **Plot Method** is **Surface**.

Addendum: Setting a specific **View Angle**, in the **View** frame, makes the system use this **View Angle** during all succeeding **Image** operations, until selecting a new **Azimuth** or **Elevation** angle again changes it. Selecting a low **Plot Density** speeds up all **Image** processing operations.

Part 10 - Single Image Window



Clicking **Single Image** in the **Analysis** menu evokes the **Single Image** window.

- 1. The **Mask** selection, **Date** and **Time** indicate how the **Image** was generated and when the **Photo** was captured. Relevant parameters are displayed in the **Statistics** panel.
- 2. Number of Cells dragging the slider to a new position, as indicated by the Edit box, selects a new number of cells used for displaying the spatial variability of the Image data points. Alternatively print a new value in the Edit box and press Enter.
- 3. **Distribution** displays the distribution of **Image Data** points using the number of cells set by the **Number of Cells** slider. The height of the bars represents the number of **Image** data points in each cell.
- 4. Cumulative displays the cumulative distribution of Image data points.
- 5. Close to close the Single Image window.

Addendum: The **Single Image** window is useful for identifying outlier **Image** data values and for the interpretation of the **Data Values** generated by the system. Multiple copies of the window can be kept open at the same time.

Part 11 - Group of Images Window



Clicking Group of Images in the Analysis menu evokes the Group of Images window.

- 1. The Mask ON/OFF, Mask selection, Subtraction mode. Mask ON/OFF, Date and Time indicate how the Image was generated and when the Photos were captured. If applicable Upper and Lower Limits are displayed.
- 2. Radio buttons in the X-Scale panel checking these radio-buttons selects the x-axis coordinates to be displayed in the diagram.
 - a. Checking **File Number** displays the file number of the file along the x-axis.
 - b. Checking Log Scale displays the values along a logarithmic x-scale.
 - c. Checking the **Image** number displays the image number of the file along the x-axis.
 - d. Checking **El. Time (Hr)** displays the elapsed time in hours, minutes and seconds along the x-axis.
 - e. Checking El. Time (Sec) displays the elapsed time in seconds along the x-axis.

The diagram is automatically updated when a radio-button is changed. The system preserves the setting of the individual check boxes until changed again by the user.

3. Check Boxes in the Statistics panel – checking these check-boxes selects the following parameters to be displayed in the diagram.

- a. Checking **SD** displays the standard deviations.
- b. Checking Max Value displays the maximum values.
- c. Checking Min Value displays the minimum values.
- 4. **Display** pull-down menu selecting diagram y-scale.
- 5. **Smoothing** pull-down menu selecting **Medium** displays all values moderately smoothed, selecting **None** displays all values un-smoothed and selecting **More** displays all values strongly smoothed.
- 6. **RGB Values** to open the **Colour Space** window and display the RGB values.
- 7. **Export Data** to export values to and Excel spreadsheet (including RGB triplet values.





- 9. **The RGBvalues Colour Space** window displays the RGB triplet values in either the RGB or the HSV colour space.
- 10. **RGB radio-button** click to display the triplets in the RGB colour space.
- 11. HSV radio-button click to display the triplets in the HSV colour space.

12. The **Start-point** displays the triplet for the first photo in the sequence. The **Stop-point** displays the triplet for the last photo in the sequence.

Part 12 - Curve Tracker Window

Clicking Curve Tracker in the Analysis menu evokes the Curve Tracker window.



- 1. If more than one **Image** is loaded, the **Curve Tracker** opens in the **Curve** mode, as displayed above.
- 2. In the **Date & Time** panel the date and time of photo capturing is displayed.
- 3. In the **White Track** panel, the statistics (**SD**, **Max Value** and **Min Value**) of the last identified region of interest (always displayed in white color), can be activated or deactivated by use of the check-boxes entitled **SD**, **Max Value** and **Min Value**.
- 4. In the **Subtract First Image** panel, checking the **Sub First Image** check-box subtracts the first image values from all the res and all curves start from the zero level.
- 5. Radio buttons in the X-Scale panel checking these radio-buttons selects the x-axis coordinates to be displayed in the diagram.
 - a. Checking **File Number** displays the file number of the file along the x-axis.

- b. Checking **Log Scale** displays the values along a logarithmic x-scale.
- c. Checking the **Image** number displays the image number of the file along the x-axis.
- d. Checking **El. Time (Hr)** displays the elapsed time in hours, minutes and seconds along the x-axis.
- e. Checking **El. Time (Sec)** displays the elapsed time in seconds along the x-axis.
- 6. When clicking the rectangular squares in the **Track Displayed** panel, the corresponding curve is temporary hidden. Clicking the actual square once more again displays the curve.
- 7. In the **Display** pull-down menu, the scale of the y-axes is set.
- 8. **Smoothing** pull-down menu selecting **Medium** displays all values moderately smoothed, selecting **None** displays all values un-smoothed and selecting **More** displays all values strongly smoothed.
- 9. Export Data button to export data to an Excel spreadsheet.
- 10. Close to close the Curve Tracker window.
- 11. If only one **Image** is loaded, the **Curve Tracker** opens in the **Statistics** mode, as displayed below.



- 12. In the **Select Mode** panel clicking the **Curves** radio-button displays the curve mode diagram.
- 13. In the Select Mode panel clicking the Bars radio-button displays the bar diagram.
- 14. Clicking the **Step Through** button displays the successive set of **ROI**s for different **Images**.

Part 13 - Gallery Window

Clicking Gallery in the Analysis menu evokes the Gallery window.

🛃 Gallery 700						_
First File: A-0001.jpg			Last File: A-	0071.jpg		
		1	22	2.2		
A-0001.jpg	A-0006.jpg	A-0011.jpg	A-0016.jpg	A-0021.jpg	A-0026.jpg	A-0031.jpg
A-0036.jpg	A-0041.jpg	A-0046.jpg	A-0051.jpg	A-0056.jpg	A-0061.jpg	A-0066.jpg
A-0071.jpg						
						Close

- 1. The **Gallery** window displays the entire sequence of **Images** in color representation as selected by **Plot Color** in the **View** frame. The layout automatically adapts to the number of **Images** included in the sequence. In total 150 **Images** can be displayed in the **Gallery** window. The actual number of the **Image** is displayed below respective **Image**.
- 2. Close to close the Gallery window.

Part 14 - Image Slicer Window

1. The **Image Slicer** window is evoked by clicking **Image Slicer** in the **Analysis** pulldown menu.



- 2. The right-most display areas (**Photo** and **Image**) in the **Profile Level** panel display the position of the actual slice.
- 3. The left-most display area displays the **Profile** of the actual slice.
- 4. The positions of the actual slice are controlled by the sliders, edit boxes and buttons within the **Horizontal** and **Vertical** panels.
- 5. Pressing the right-most or the left-most arrow in the slider moves the slice to the right (Vertical) or up (Horizontal) respectively to the left (Vertical) or down (Horizontal). The actual position is displayed in the Edit box.
- 6. Dragging the slider moves the slice directly to the target position.
- 7. Clicking **Scan** makes the slice scan through the complete image.
- 8. The **Scan** procedure can be halted at any time by clicking the **Pause** button. The **Scan** procedure is continued by clicking the same button now labeled **Cont**.
- 9. The **Scan** procedure can be terminated at any time by clicking the **Stop** button.

- 10. The **Display** pull-down menu can be set to either **Fix**, **Auto** or **Fit**. **Fix** mode implies that the Y-axis is fixed with the end-points set to the maximum and minimum value respectively. **Auto** mode implies that the upper end-point of the Y-axis adapts to the maximum value of the actual **Profile**. **Fit** mode implies that both the upper and lower end-point of the Y-axis adapt to the actual **Profile** maximum and minimum values.
- 11. Clicking the **Hold** button implies that each new **Profile** is superimposed on the previously displayed **Profiles**. Clicking the same button again, now labeled **Erase**, returns the system to displaying only a single **Profile**.
- 12. Clicking the **Close** button closes the **Image Slicer** window.
- 13. The **File Name**, **Date** and the **Time** of **Image** capturing are displayed in the left side of the **Control Panel**.
- 14. The **Image** size as well as the **Maximum** and **Minimum** value of the actual **Profile** are displayed in the top right corner of the **Profile Level** window.

Part 15 - Cross Sectioner Window

To display a sequence of cross sections, a specific **Cross Section** of the **Images** needs to be loaded. Proceed as follows.

- 1. Load a **Photo** sequence (in this case the test image sequence **A-0001** to **A-0071** stored in the **MN DATA** folder).
- 2. Select **Cross Sectioner** from the **Analysis** pull-down menu. The **Cross Sectioner Prepare 700** window will now open.



- 3. Select Vertical Line or Horizontal Line by checking the associate radio-button.
- 4. Click somewhere in the **Image** to generate the line.
- 5. Click Accept to accept the line selection. The Cross Sectioner 700 window will no open.



- 6. Date & Time panel displays date and time of photo capturing.
- 7. **X-Scale** panel check radio-button to select t x-scale.
 - a. Checking **File Number** displays the file number of the file along the x-axis.
 - b. Checking the **Image** number displays the image number of the file along the x-axis.
 - c. Checking **El. Time (Hr)** displays the elapsed time in hours, minutes and seconds along the x-axis.
 - d. Checking El. Time (Sec) displays the elapsed time in seconds along the x-axis.
- 8. Display panel check the Hide Grid and Hide Axis to hide the grid and/or the axis.
- 9. Control Panel
 - a. Check **Top View**, **Side View**, **Perps 1**, **Persp 2**, **Persp3** or **Persp 4** to select the view angle.
 - b. Drag the **Color** slider to alter the diagram color scale.
 - c. Select **Plot Method** in the **Plot Method** pull-down display.
 - d. **Export Data** to export data to Excel.

- e. **Close** to close the window.
- **10. First** and **Last File Name, Cross Section Number** and **Mode** are displayed above the **Cross Section Map** panel.

Part 16 - Profiler Window

The **Profiler** is useful for displaying and analyzing circular erythema and blanching reactions.

- 1. Load a **Photo** sequence (in this case the test image sequence **UVB-0001** to **UVB-0004** stored in the **PROFILER FILES**.
- 2. Select **Profiler** form the **Analysis** pull-down menu. Click the **Show Last Photo** radiobutton. The **Profiler** window should now look like:



- 3. Click the **Place High** button. After a while 5 concentric rings will be automatically placed around the peak value of the reaction.
- 4. Drag the **Number of Rings** slider to the right until the edit box displays 10. Drag the **Size of Rings** slider to the right until the edit box displays 130. The **Profiler** window should now look like (alternatively insert the number directly in the edit box and press **Enter**) :



5. Select **Curve** under the **Plot Method** popup menu. The **Profiler** window should now look like:

TiVi700 Analyzer Manual



- 6. The **Curve** displays the average TiVi value within the different rings and the associated standard deviation.
- 7. Click the **Step Through** button repeatedly to step through the sequence of photos uploaded.
- 8. Click the **Display Curves** button to display all curves in a separate window:



- 9. In the **Control** panel:
 - a. Show SD to display or hide the standard deviation bars.
 - b. **Display** to set the y-axis.
 - c. **Display Bars** to display the result as a bar graph.
 - d. **Close** to close the window.
- 10. Click the **Display Bars** button to display the result as a bar-graph.
 - a. Display Panel
 - i. Mean Intensity to display the mean intensity.
 - ii. Mean Intensity & SD to display standard deviations on the mean value bars.
 - b. Mode Panel
 - i. Curve Number displays curve number along x-axis.
 - ii. Ring Number ring number along the x-axis.
 - c. Bar Graph Panel
 - i. 2D displays the diagram in 2D mode.
 - **ii. 3D** displays the diagram in 3D mode (can be rotated).
 - d. Pallette pull-down menu -to select bars color palette.
 - e. Close button to close the window.
- 11. **Export Data –** exports data to an Excel spreadsheet.



Part 17 - Wizard

The **Wizard** is intended for rapid and automatic analysis of a sequence of **Images** in which **ROIs** have been drawn. The **Wizard** calculates the average TiVi-value inside the **ROIs** including only pixels above or below set threshold.

First the **Welcome to the Express Wizard** is displayed in which the **Above Threshold** or **Below Threshold** mode can be selected and the **Threshold** level set. When clicking the OK button these values are accepted and the **TiVi Chart** window is displayed

The TiVi Chart window displays the calculated TiVi values for each ROI and Image.



With Threshold = 50% and **Above** mode selected.



With Threshold = 50% and **Below** mode selected.

Part 18 – TiVi700 Chart Window

- 1. Mean Values radio-button check to display average TiVi-values.
- 2. Median Values radio-button check to display median TiVi-values.
- 3. Area radio-button check to display the surface area of individual patches.
- 4. **Mean * Area** check to display the product of average TiVi-Values and the surface area.
- 5. **Sub First Image check-box** check to subtract the first Image **ROI** TiVi-values from the remaining TiVi-values in the sequence of **Images.**
- 6. **Standard Dev check-box** check to display the standard deviations of the last (white) ROIs.
- 7. Average Curve check to display the average curve.

- 8. Average Std. check to display the standard deviation of the average curve.
- 9. Select Reference Curve subtracts the TiVi-values in the selected curve from al the remaining TiVi-values. This feature is useful to relate TiVi-values of various **ROIs** to those of a control **ROI**.
- 10. First File No. Displayed to select the first file number displayed.
- 11. Last File No. Displayed to select the last file number displayed.
- 12. **Smoothing** to select the degree of smoothing.
- 13. Select X-axis to select the x-axis scale.
- 14. **Select Y**-axis to select the y axis scale.
- 15. % of Max Value creates vertical lines at Images when the set percentage of maximum value (indicated by a dot) has been reached.
- 16. Bar Graph to display the result in terms of a bar graph.



- a. Display panel.
 - i. Mean Int displays the average TiVi values.
 - ii. **Mean Int & SD** displays the average TiVi-values with standard deviations.
 - iii. Median Int displays the median TiVi-values.
 - iv. Area displays the surface area of the patches

- b. Mode panel
 - i. **File No. –** sets the file number along the x-axis.
 - ii. **ROI** sets the ROI number along the x-axis.
- c. Bar Graph panel
 - i. **2D** to display a two-dimensional bar graph.
 - ii. 3D to display a three dimensional bar graph. By placing the mouse cursor inside and pressing the left mouse button, the BarGraph panel the 3D view can be altered.
- d. **Pallete** to select the color scale of the bar graph
- e. **Close** to close the ChartBG700 window.
- 17. Close to close the TiVi700 Chart window

18. ROIs Displayed panel

- a. **Color buttons –** click a color button to temporarily hide the associated curve.
- b. Edit box displays the time/file number/image number at which vertical lines crosses the curves.

19. File Menu

- a. Save TiVi Chart to save the actual TiVi700 chart (can be re-opened from the TiVi700 Analyzer window.
- b. **Close –** to close the actual TiVi Chart.
- c. Export to export data from the TiVi Chart to an Excel spreadsheet.

Above the panels information about the **First File Name, Last File Name, Step, Number of Files, Threshold setting** and **Mode** is displayed.

Part 19 - Project Window

The **TiVi700 Project** window displays average TiVi-values by importing these values from a saved **TiVi700 Chart**. While the **TiVi700 Chart** generally represents data recorded from a single subject, the **TiVi700 Project** displays average data values from up to 20 subjects.

A new **TiVi700 Project** window is opened from the **TiVi700 Analyzer** window (File→Project→ New TiVi Project).

An existing **TiVi700 Project** window is opened form the **TiVi700 Analyzer** window (File \rightarrow Project \rightarrow Open TiVi Project).

In the **TiVi Project** window below, seven average TiVi value curves have been uploaded from saved TiVi Charts (AA_CHART \rightarrow GG_CHART).



- 1. Objects panel lists all uploaded objects.
 - a. Remove object permanently by clicking on its corresponding quadrat in the **Objects** panel
- 2. Controls panels
 - **a.** Mean radio-button click to display average values.
 - b. Median radio-button click to display median values.
 - c. Area click to display the average patch area values.
 - d. Mean * Area Click to display the product of Mean an surface area values.
 - e. Average check box check to display the average of all curves displayed.
 - **f.** Standard Deviation check box check to display the associated standard deviations.
 - g. Log Scale check-box click to display data in logarithmic scale.
 - h. Smoothing to smooth the curves displayed (Medium, None, More).
 - i. Select Y-axis to set scale of Y-axis.
 - **j.** % of Max Value to display vertical lines corresponding to the file number at which the set percentage of the maximum value (indicated by a square) is reached. The file number is also displayed in the color quadrat next to the name in the **Objects** panel.
- 3. File Menu

- a. Save TiVi Project to save the actual TiVi Project.
- b. Export Data to export data displayed to an Excel spreadsheet.
- c. Close to close the TiVi700 Project window.
- 4. Load Objects to upload a new object to the TiVi700 Project window. Click on a saved TiVi700 Chart to upload the data.
- 5. Close to close the TiVi700 Project window.

Part 20 - Excel Files Generated

The following Excel Files are generated by the various windows of the TiVi700 Analyzer:

TiVi700 Analyzer	Main Window:
Export Image:	Full Image.
Group of Images Export Data:	Window: Mean Intensity Values, Standard Deviations, Min Intensity Values, Max Intensity Values
Curve Tracker W Export Data:	indow: Mean Intensity Values, Median Intensity Values, Standard Deviations, Min Intensity Values, Max Intensity Values
Image Slicer Wind	dow:
Export Data:	Actual profile.
Cross Sectioner W	Vindow:
Export Data:	Full Image
Profiler Window (Display Curves):
Export Data:	Mean Values, Standard Deviations, Points
Chart Window: Export Data:	Mean TiVi Values, Extension.
Project Window:	Mean TiVi Values, Standard Deviations, Extension.
Export Data:	Perc-F_T: percentage, Image No for % of Max for individual ROIs

Part 21 – Specifications

Property	Value
Туре	
Camera Head	Canon EOS 550D / 650D/ 700D
Computer	64 bit, 4GB min (8GB recom.)
Flexible Arm	Neo-Flex
Dimensions	
Computer	4x26x36 cm, 3 kg
Camera	5x18x12 cm, 0.5 kg/ 24x16x13, 1 .5 kg
Flexible Arm	about 50 cm, 1.2 kg
Photo format	jpeg
Photo size (pixels)	1728 x 2592 - 3456 x 5184
Maximal resolution (at 30 cm)	50 micrometers
Distance	10 cm – inf.
Field of View (at 30 cm)	280 x 200 mm
Maximal photos per minute	12
Maximal number of photos in sequence	9999
Instability	
Temperature drift	0.35% per Degrees C
Drift over 2 months period	< 1%
Image uniformity in sensitivity	
Without normalization	3% (corners) 2% (centre)
With normalization	0.25% /corners) 0.13% (centre)
Distance dependence	Not measurable $(13 - 28 \text{ cm})$
Ambient light dependence	
40W bulb at 50 cm distance	10%
11W fluorescent light source at 50 cm distance	10%

Part 22 - Demo Assistant

The *Demo Assistant* is an online interactive guide to the effective use of the basic functions and features of the TiVi700-system. The *Demo Assistant* is available from the **Demos Assistant** pull-down menu. When a topic in the *Demo Assistant* is selected a specific topic demo window is opened which guides the user step-by step on how to perform certain tasks.

TiVi700 Analyzer 1.1.1					Welcows to the Londing Dhat
File Edit Page Analysis Wizard Maintenar	ice Help Language	Tool Boxes Demo Assista	nt	r	welcome to the toading Phot
🖬 🕫 🖬 🖬 🗭 🖬 🔺 💻 🖬 🕿 🕿	Mask OFF	lanan	S	how / Crop	This demo shows how to load a sequence of photos with a user defined step, change the color scale and step through the sequence of photos automatically and
ROIs: PetCR Dist(mm).0 Zoom0	Photo Size: 75375 Date: 21-Feb-2008 Time First: 09:20:36 Time: 09:21:40 00:01:04 Time Last: 09:27:40	LowerLimt0	Upper Limit 100 T/V	Mean 500 119 450 115 400 SEM 400 0 350 SD 1 300 Mn 44 250 Max 200 223/223 Width 150 225 100 Height 335 50 Points 0 75553	manually. Click the First Photo button and navigate to the MI DATA folder located in in the TW/600demonstration folder. Double-click the A-0001 file. Click the Last Photo button. Double-click the A-0061 file. Click the Step button several times until the Select Step edit box displays 10. Set Color Span to 500. Click the Automatic button to scan through the photos.
Actual Photo Actual Photo A-0011.jpg	Show First Photo	Reverse	Pause	Comments	Click the Manual button twice to display the first two photos in the uploaded sequence.
First Photo A-0001.jpg	Show Last Photo First Photo	Manual Reset	Automatic Abort	Color Span 500 Delay (sec) 0	This completes the demo. Click the Close button to close this window.
Last Photo A-0071.jpg Subtracted Photo	Last Photo	Crop Images-	Ma Subtract Fifr	ask x1: x2 ↓: y1: y2:	
Clear All Reverse	10 Step	Overlay Crop & Save	Upper Limit	Reset Mask Refresh	Ciose

The Demo Assistant includes to following topics:

DEMO ASSISTANT TiVi File Names Loading Photos Crop & Save Thresholding Movies **Image Subtraction ROI** Management Draw ROIs Delete ROIs Move ROIs Stretch ROIs Set ROI Parameters Save & Load ROIs Analysis Single Image Groups of Images

Curve Tracker Gallery Image Slicer Cross Sectioner Profiler Wizard TiVi Chart Project TiVi Pre-Processor Using the *Demo Assistant* is an effective way of getting acquainted with the basic functions and features of the TiVi700 system.

Part 23 DropData

The **Dropdata** feature is a convenient way of bringing image data in focus directly into a spread-sheet by a single click on a spread-sheet cell. The associated **Organizer** holds editable subject data and the subject identifiers can be exported directly to the **DropData** spread-sheet. Both the **DropData** spread-sheet and the **Organizer** spread-sheet are fully editable and both can be saved as **Excel** documents for further processing.

The DropData / Organizer features are best explained by an example.

Assume that the response to topical application of methyl nicotinate from three separate subjects has been stored in the files **A-0001 – A-00071**, **C-0010 – C-0085** and **E-0010 – E-0095** stored in the *MN DATA* folder in the *TiVi700 demonstration* folder.

- 1. Click the **First Photo Photo** button in the **TiVi700 Analyzer** window and select the A-0001 file in the *MN DATA* folder located in the *TiVi700demonstration* folder.
- 2. Select **Data to Sheet** under the **DropData** pull-down menu. The **DropDataSettings** window will now open.

DropDataSettings	_ ×
Number of Rows	20
Number of Columns	10
	ОК

3. Click the **OK** button to open the **DropDataSheet700** window.

TiVi700 Analyzer Manual

					Aean Valu	es				
		Upper L	.imit: 100	Mask	: FiFr	Nar	ne			
		Lower	Limit: 0	Whole	e Image	Pho	oto: A-007	1		
	A	В	С	D	E	F	G	н	I	J
1	0	0	0	0	0	0	0	0	0	
2	0	0	0	0	0	0	0	0	0	
3	0	0	0	0	0	0	0	0	0	
4	0	0	0	0	0	0	0	0	0	
5	0	0	0	0	0	0	0	0	0	
6	0	0	0	0	0	0	0	0	0	
7	0	0	0	0	0	0	0	0	0	
8	0	0	0	0	0	0	0	0	0	
9	0	0	0	0	0	0	0	0	0	
10	0	0	0	0	0	0	0	0	0	
11	0	0	0	0	0	0	0	0	0	
12	0	0	0	0	0	0	0	0	0	
13	0	0	0	0	0	0	0	0	0	
14	0	0	0	0	0	0	0	0	0	
15	0	0	0	0	0	0	0	0	0	
16	0	0	0	0	0	0	0	0	0	
1/	0	0	0	U	0	0	0	0	U	
18	0	0	0	0	0	0	0	0	0	
19	0	0	0	0	0	0	0	0	0	
20	0	0	0	0	0	0	0	0	0	

4. Move the **DropDataSheet700** window to a location next to the **TiVi700Analyzer** window.

TiVi700 Analyzer 1.1.1 To be used with Camera wit	th Serial Number: 22333	50670		×	٦ŗ.	DropDat	aSheet700							-	□ X
File Edit Page Analysis Wizard DropData I	Maintenance Help L	anguage Tool Boxes Demo Assi	tant		2										
9 8 8 9 9 1 4 4 4 2 2 2	 3 					- Tabl	e								
A19. 7D .	, <u> </u>		Show /	Cron						Me	an Values				
Whenle Bridge ROI OFF Subtract OFF	Mask OFF		0.1017	orop				Upper L	.imit: 100	Mask: I	FIEr	Name:			
Photo	Dhata Sizar	Image		Have	1.										
Rois. Polick Dist(min).0 200m.0	75375	Lower Link.0 Opper Link	.100 1101	- 700 223				Lower I	Limit: 0	Whole I	mage	Photo: A-00	71		
Sector Contraction of Sector				Median											
				630 225			A	B	C	D	E F	G	н	I	1
				560 SEM		1	0	0	0	0	0	0 0	0	0	
				0		2	0	0	0	0	0	0 0	0	0	
ALC: NOT THE REPORT OF THE REPORT OF				490 35		4	0	0	0	0	0	0 0	0	0	0
and the second				420 Min		5	0	0	0	0	0	0 0	0	0	0
	Date:			350 0/0		6	0	0	0	0	0	0 0	0	0	0
	21-Feb-2006			Max		7	0	0	0	0	0	0 0	0	0	0
	Time First:		Sec. and	280 521/521		8	0	0	0	0	0	0 0	0	0	
	Time:			210 225		10	0	0	0	0	0	0 0	0	0	
	09:27:40			Height		11	0	0	0	0	0	0 0	0	0	0
	00:07:04			140 335		12	0	0	0	0	0	0 0	0	0	0
	1 me Last: 09:27:40			70 Points		13	0	0	0	0	0	0 0	0	0	0
and the second				0 75358		14	0	0	0	0	0	0 0	0	0	0
Palast Photon		Stop Through			1	15	0	0	0	0	0	0 0	0	0	
		Step mough				10	0	0	0	0	0	0 0	0	0	
Actual Photo	Show First Photo	Reverse	Pause	Comments		18	0	0	0	0	0	0 0	0	0	0
A-0071.jpg	Show Last Photo	Manual	Automatia			19	0	0	0	0	0	0 0	0	0	0
First Photo			Automatic	Color Span 700		20	0	0	0	0	0	0 0	0	0	0
A-0001.jpg	First Photo	Reset	Abort	Delay (sec)											
Last Photo		One la contra de l			1									-	
A-0071.jpg	Last Photo	Crop images	Mas	sk	11	[Cont	trois—							_	
Subtracted Photo	Step	Show ROI S	ibtract FFr	x xt x2 y yy yz		Clic	k on a Cell in t	the Table	to insert d	ata from Ima	ge				
	5	Up	er Limit	Denet Mark			Load Table	0	Show M	eans	Set Row Na	mes Ba	Graph	Res	et
			100	Reset Mask		-									
Clear All Reverse	Step	Crop & Save	ver Limt	Refresh		1	Save Table	6	Show Si	D	Set Col Nan	nes Exp	ort Data	Clos	e
					1			-							

- 5. Click the **Show First Photo** radio-button in the **TiVi700 Analyzer** window to display the first photo of the sequence.
- 6. Click cell A1 in the **Dropdata700** spread-sheet to insert the *Image Mean value* in focus into the **DropData700** table.
- 7. Click the **Manual**-button in the **TiVi700 Analyzer** window to display the next photo / image in the sequence. Click cell B1 in the **DropData700** spread-sheet to insert the *Image Mean value* in focus into the DropData700 table.
- 8. Click the **Show Last Photo** radio-button to display the last photo of the sequence.
- 9. Click the **Reverse**-button and then the **Manual** button and then again click the **Forward** button to display photo A-0066.
- 10. Click cell C1 in the **Dropdata700** spread-sheet to insert the *Image Mean value* in focus into the **DropData700** table.
- 11. Click the **Manual**-button in the **TiVi700 Analyzer** window to display the next photo / image in the sequence. Click cell D1 in the **DropData700** spread-sheet to insert the *Image Mean value* in focus into the DropData700 table.
- 12. The DropData700 spread-sheet should now look like:

TiVi700 Analyzer Manual

				1	Aean Valu	ies				
		Upper L	.imit: 100	Mask	: FiFr	Nar	me			
		Lower	Limit: 0	Whol	e Image	Pho	oto: A-006	51		
	A	В	С	D	E	F	G	н	I	J
1	119	119	206	215	0	0	0	0	0	
2	0	0	0	0	0	0	0	0	0	
3	0	0	0	0	0	0	0	0	0	
4	0	0	0	0	0	0	0	0	0	
5	0	0	0	0	0	0	0	0	0	
6	0	0	0	0	0	0	0	0	0	
7	0	0	0	0	0	0	0	0	0	
8	0	0	0	0	0	0	0	0	0	
9	0	0	0	0	0	0	0	0	0	
10	0	0	0	0	0	0	0	0	0	
11	0	0	0	0	0	0	0	0	0	
12	0	0	0	0	0	0	0	0	0	
13	0	0	0	0	0	0	0	0	0	
14	0	0	0	0	0	0	0	0	0	
15	0	0	0	0	0	0	0	0	0	
16	0	0	0	0	0	0	0	0	0	
17	0	0	0	0	0	0	0	0	0	
18	0	0	0	0	0	0	0	0	0	
19	0	0	0	0	0	0	0	0	0	
20	0	0	0	0	0	0	0	0	0	

- 13. The mean TiVi-values of the two first and two last images are now inserted and displayed in the **DropData700** spread-sheet.
- 14. Repeat the procedure with photo sequences C and E.
- 15. The **DropData700** spread-sheet should now look like:

	Upper L	imit: 100			2220				
			Mask	:: FiFr	Nai	ne			
	Lower I	.imit: 0	Whol	e Image	Ph	oto: E-009	5		
A	В	C	D	E	F	G	н	I	J
1 119	119	206	223	0	0	0	0	0	(
2 153	155	201	217	0	0	0	0	0	
3 142	136	185	186	0	0	0	0	0	(
4 0	0	0	0	0	0	0	0	0	(
5 0	0	0	0	0	0	0	0	0	1
6 0	0	0	0	0	0	0	0	0	3
7 0	0	0	0	0	0	0	0	0	
8 0	0	0	0	0	0	0	0	0	
9 0	0	0	0	0	0	0	0	0	1
10 0	0	0	0	0	0	0	0	0	1
11 0	0	0	0	0	0	0	0	0	
12 0	0	0	0	0	0	0	0	0	
13 0	0	0	0	0	0	0	0	0	
14 0	0	0	0	0	0	0	0	0	(
15 0	0	0	0	0	0	0	0	0	(
16 0	0	0	0	0	0	0	0	0	
17 0	0	0	0	0	0	0	0	0	0
18 0	0	0	0	0	0	0	0	0	1
19 0	0	0	0	0	0	0	0	0	(
Second	0	0	0	0	0	0	0	0	(

TiVi700 Analyzer Manual

				I	Aean Valu	es				
		Upper L	.imit: 100	Mask	:: FiFr	Nai	ne			
		Lower	Limit: 0	Whol	e Image	Ph	oto: E-009	5		
	A	В	С	D	E	F	G	н	I	J
1	119	119	206	215	0	0	0	0	0	
2	153	155	153	155	0	0	0	0	0	
3	142	136	185	186	0	0	0	0	0	
4	0	0	0	0	0	0	0	0	0	
5	0	0	0	0	0	0	0	0	0	
6	0	0	0	0	0	0	0	0	0	
7	0	0	0	0	0	0	0	0	0	
8	0	0	0	0	0	0	0	0	0	
9	0	0	0	0	0	0	0	0	0	
10	0	0	0	0	0	0	0	0	0	
11	0	0	0	0	0	0	0	0	0	
12	0	0	0	0	0	0	0	0	0	
13	0	0	0	0	0	0	0	0	0	
14	0	0	0	0	0	0	0	0	0	
15	0	0	0	0	0	0	0	0	0	
16	0	0	0	0	0	0	0	0	0	
17	0	0	0	0	0	0	0	0	0	
18	0	0	0	0	0	0	0	0	0	
19	0	0	0	0	0	0	0	0	0	
20	0	0	0	0	0	0	0	0	0	

- 16. Click the **Show SD** radio-button to display the standard deviation values in the **Table**.
- 17. Click the **Clear All** button to clear all edit boxes.
- 18. Click the **Set Row Names** button and insert the names of the files in the corresponding edit boxes.

📣 SetRov	vs700	- • ×
Raw 1	A	Raw 2 C
Raw 3	E	Raw 4
Raw 5		Raw 6
Raw 7		Raw 8
Raw 9		Raw 10
Raw 11		Raw 12
Raw 13		Raw 14
Raw 15		Raw 16
Raw 17		Raw 18
Raw 19		Raw 20
		ок

19. Click the OK button. The DropDataSheet 700 should now look like:

			Mean Values							
		Upper l	.imit: 100	Mask	: FiFr	Nai	me:			
		Lower	Limit: 84	Whole	e Image	Ph	oto: E-009	5		
	A	В	с	D	E	F	G	н	I	J
Α	119	119	206	223	0	0	0	0	0	
С	153	155	201	220	0	0	0	0	0	0
E	142	136	186	186	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	
	0	0	0	0	0	0	0	0	0	
	0	0	0	0	0	0	0	0	0	
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	
	0	0	0	0	0	0	0	0	0	(
	0	0	0	0	0	0	0	0	0	

20. Click the **Set Col Names** button, click the **Clear All** button and insert the appropriate column names. Click the **OK** button.

Upper	Limit: 100							
		Mask:	FiFr	Name	•			
Lower	Lower Limit: 0 Whole Image				Photo: E-0095			
first2	last1	last2						
119	206	223	0	0	0	0	0	3
155	201	217	0	0	0	0	0	
136	185	186	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
			0	0	0			
	first2 119 155 136 0	first2 last1 119 206 155 201 136 185 0 0	Iost1 Iost2 Iost2 119 206 223 155 201 217 136 186 137 185 186 0 0 0	first2 last1 last2 and 119 206 223 00 155 201 217 0 136 185 186 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	first2 last1 last2 119 206 223 0 0 155 201 217 0 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 <td< td=""><td>first2 last1 last2 119 206 223 00 00 0 155 201 217 0 0 0 136 185 186 0 0 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0</td><td>first2 last1 last2 o o o o 119 206 223 0 0 0 0 155 201 217 0 0 0 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 <</td><td>first2 last1 last2 o o o o 119 206 223 0 0 0 0 0 0 155 201 217 0 0 0 0 0 0 136 188 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 <td< td=""></td<></td></td<>	first2 last1 last2 119 206 223 00 00 0 155 201 217 0 0 0 136 185 186 0 0 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	first2 last1 last2 o o o o 119 206 223 0 0 0 0 155 201 217 0 0 0 0 136 185 186 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 <	first2 last1 last2 o o o o 119 206 223 0 0 0 0 0 0 155 201 217 0 0 0 0 0 0 136 188 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 <td< td=""></td<>

- 21. Click the **Save As** button and save the table under a selected name and folder.
- 22. Mark the cells holding numerical values and click the **Reset** button to empty the spread-sheet of data.
- 23. Click the Load Table button and load the saved table again.
- 24. Click the **Bar Graph** to display data as a Bar Graph.



- 25. The graph can be rotate by placing the mouse pointer in the diagram, pressing the left mouse button and moving the mouse. Click the **Close** button.
- 26. Click the **Export Data** button to save the Table as an **Excel** document.
- 27. Click the Close button to close the DropDataSheet700 window.

Detailed description of **DropData**:

DropData settings:

- 1. Number of Rows edit box to insert the actual number of rows in the DropData Sheet700 (up to 20).
- 2. Number of Columns edit box to insert the actual number of columns in the DropDataSheet700 (up to 10).
- 3. OK button to open a new DropDataSheet700

DropDataSheet700:

 column heading, pressing the left mouse button and moving the point to a new position.

- 2. Values above the spread-sheet displays the actual image **Upper Limit** and **Lower Limit** settings, the **Mask** setting and if the image is cropped as wells as the name of the actual photo and the name of the **DropDataSheet700** (if saved).
- 3. Save As button to save the Table under a user selected name.
- 4. Save button to save the Table under the name selected with Save As.
- 5. Load Table button load a saved Table.
- 6. Show Means radio-button to display the actual image mean values.
- 7. Show SD radio-button to display the actual image standard deviations.
- 8. Set Row Names to open the Row editor window in which names of the rows can be inserted.
- 9. Set Col Names to open the Column editor window in which names of the columns can be inserted.
- 10. Bar Graph to open the Bar Graph window.
- 11. Export Data to export spread-sheet data to an Excel document.
- 12. Reset to reset the content of the table. Select an area by pointing at the first cell.
 Press the mouse button and the drag the mouse to a new cell thereby defining the area.
 Release the mouse button. The defined area of cells now is displayed in blue background colour. Click the Reset button. Alternatively cells can be selected by use of the Shift tangent. If the content of a single cell is to be deleted, point at the cell and print "0" in the cell edit box.
- 13. Close to close the DropDataSheet700.

Part 24 Organizer

1. Select **Organizer** under the **DropData** pull-down menu. The **Organizer 700** welcome window open.

📣 Orga	anizer 700	x
	Welcome to The Organizer	
	New Organizer	
	Open Organizer	
	Close	

2. Click the **New Organizer** button to open a new **Organizer**.

New Organizer 700 💻 😐 🗙	
Substance Tested	
Comments	
•	
User Defined Fields	
Create New and Save	
Close	

- 3. Print *Methyl Nicotinate* in the **Substance Selected** edit box.
- 4. Print *This is my first Organizer project* in the **Comment** area
- 5. Click the User Defined Fields button to open the User Defined Fields window.

📣 Ueser Defined Fi 💻 🗖 🗙
Field 1
Field 2
Field 3
Just Close Save & Close

6. Print 'Additional Comments' in the Field 1 edit box. Leave the other edit boxes blank.

Ueser Defined Fi 💻 😐 🗙
Field 1
Additional Comments
Field 2
Field 3
Just Close Save & Close

- 7. Click the Save & Close button.
- 8. In the New Organizer 700 window click the Create New and Save button.

9. Save the **Organizer** under a selected name (e.g. T4) and folder. **The Organizer 700** named T4 now opens.

ime: T4	Substance Tested: Methyl Nicotinate									
ID	Subject	Gender	Date of Birth	Address	Phone No.	e-mail	Test Date	Tester	Done	Additional Comme
<										

10. Click the **Add New Subject** button to add data for a new subject to the T4-table. The **Add New Subject 700** now opens. Fill out Subject data.

🛃 Add N	lew Subject 700		- 🗆 X
- Sut	biect Data		
Sul	bject Identifier	MA	
Ge	nder	Male Fem	ale
Dat	te of Birth	1948-02-23	
Sul	bject Address	Main Street 100	
Sul	bject Phone	01203-235456	
e-n	nail	martin@glocalnet.net	
Tes	sting Date & Time	2012-02-25-12-20	
Tes	ster Identifier	Theresa M	
Ad	ditional Comments	Not tested befow	
	Save & Close	Just C	lose

11. Click the Save & Close button. The Organizer 700 window should now look like:

Additional Comm Not tested befow
,

12. Fill out data for two more subjects accordingly.

Organize	er 700										
Name: T4		Substance Tested: Methyl Nicotinate									
	ID	Subject	Gender	Date of Birth	Address	Phone No.	e-mail	Test Date	Tester	Done	Additional Comme
	1	MA	М	1948-02-23	Main Street 1	01203-235456	martin@gloc	2012-02-25	Theresa M		Not tested befow
	2	SB	F	1963-06-13	Lower Street	0203-243561	sally@hotmai	. 201202-25-1	Theresa M		
	3	CC	M	1961-04-21	Upper Street	0203-294356	charlie@tele	2012-12-25	Samuel B		excema
4						m					•
•		Remo	ove Subject			m					4

- 13. Click the **Save Organizer** to save the **Organizer** spread-sheet.
- 14. Click the **Export Data** button to export data to an **Excel** spread-sheet.

- 15. Click the **Close** button to close the **Organizer** window.
- 16. In the **TiVi700 Analyzer** window, click the **First Photo** button and load the *A*-0001.jpg file.
- 17. Select **Organizer** from the **DropData** pull-down menu.
- 18. In the **Organizer700** window, click the Open Organizer and select the Organizer just saved (T4).

ID Subject 1 MA 2 SB	Gender Date o M 1948-02	of Birth Address	Phone No. e-mail			1
1 MA 2 SB	M 1948-0			lest Date les	ter Done	Additional Comme
2 SB		02-23 Main Street 1	01203-235456 martin@gloc	2012-02-25 Theres	a M 📃	Not tested befow
	F 1963-04	06-13 Lower Street	. 0203-243561 sally@hotma	i 201202-25-1 Theres	a M 📃	
3 CC	M 1961-0	04-21 Upper Street	. 0203-294356 charlie@tele	2012-12-25 Samuel	B	excema

19. Click the **IDs to Sheet** button to open the **DropData 700** window with the subject identifiers imported.

	Mean Values														
	UL	UL			MF				Name						
	ш	LL		ShowCrop				filename			B				
		A	В	С		D	E		F	G		н	I	J	
	MA	0	0		0	0		0	0		0	0	0		0
	CC	0	0		0	0		0	0		0	0			0
Controls	Cell in the	Table t	o insert	data f	rom li	mage									

- 20. Place the DropData700 window next to the TiVi700 Analyzer window.
- 21. Click Set Col Names button and insert appropriate column names.
- 22. Transfer image data from the **TiVi700 Analyzer** to the spreadsheet (see section **DropData** above).

Detailed description of **the Organizer**:

Organizer700 settings window:

- 1. New Organizer to open a new Organizer.
- 2. Open Organizer to open a saved Organizer.
- 3. Close to close the Organizer700 window.

New Organizer 700 window.

1. Substance Tested - insert the name of the substance tested here.

- 2. **Comment** insert comments regarding the study.
- 3. User Defined Fields in addition to the standard fields in the Organizer, three optional fields can be inserted. The names of these fields are inserted in the edit boxes of the User Defined Fields window.
- 4. Create New and Save to create and save a new Organizer.
- 5. Close to close the New Organizer 700 window.

Organizer 700:

- 1. **Table** holds data regarding subject ID (number) and Subject Identifier, Gender, Date of Birth, Address, Phone Number, e-mail, Test Data, Tester ID, Done (to check when the test for the actual subject is complete) and up to three user-defined fields. The table is fully editable.
- Add New Subject opens window for including a new subject in the Organizer. This is done through the Add New Subject 700 window in which data of Subject Identifier, Gender, Date of Birth, Subject Address, Subject Phone Number, e-mail, Testing Date and Time, Tester Identifier n additional user-defined fields. Clicking Save & Close automatically insert the data into the Organizer.
- 3. **Remove Subject** to remove a subject from the **Organizer**. The unique subject ID number is not changed by this procedure.
- 4. Export Data to export the content of the Organizer to an Excel spread-sheet.
- 5. **IDs to Sheet** to export the subject ID data to a new **DropDataSheet** which is automatically generated and opened.
- 6. Edit Comments to display and edit comments related to the project.
- 7. Save Organizer saves the current content of the Organizer.
- 8. Close to close the Organizer.

Part 25 Software Installation Procedure

If another version of TiVi700Analyzer is already installed, this must first be de-installed. Go to the Control Panel \rightarrow Program (de-install a program) and select TiVi-1.1 Tissue Viability Imager (about 1GB) (or TiVi700Analyzer). Check "Yes" to de-install.

Move the folder named **TiVi700Analyzer** to the install computer desktop from the installation USB memory stick.

- 1. Click the **TiVi700 Analyzer** software setup icon in the **TiVi700Analyzer** folder. The installer will start to extract the **TiVi700Analyzer.msi** file.
- 2. Click the Next button in the Welcome to the InstallShield Wizard for TiVi700Analyzer.
- **3.** Read the *License Agreement*. Check the "*I accept the terms in the License agreement*" radiobutton to continue with the installation procedure.
- 4. Click the Next button.
- **5.** Click the Install button to install the software in the C:\Program Files (x86)\Wheelsbridge\TiVi700Analyzer folder.
- 6. When the User Control Icon is displayed click on it.
- 7. In the *User Account Control* click "*Yes*". The files will start to extract (will take about a minute or two.
- **8.** In the TiViAnalyzer InstallShield Wizard check the Launch the Program checkbox (if not already checked. Click the Finish button.
- **9.** IN the User Account Control click the "Yes" button to start installation of the MATLAB libraries. The Winzip Self-Extractor MCRINstaller.exe starts to extract the library files.
- **10.** In the *MATLAB Runtime Installer* click the **Next** button.
- **11.** Click the "*Yes*" radio button in the *License Agreement* and the click the **Next** button to start the installation of the MATLAB libraries.
- **12.** Click the Next button in the Folder selection window.
- **13.** Click "Yes" in the Folder Selection window as answer to the question " Would you like to create it?"
- **14.** Click the **Install** button in the *Confirmation* window to start the installation (will take several minutes).
- 15. Click the Finish button to complete the installation process.
- 16. To start the TiVi700 Analyzer click the Launch TiVi700Analyzer icon on the desktop.

If the program does not start correctly the *Wheelsbridge* window is probably not set up for reading and writing by the User. Proceed as follows to set the *Wheelsbridge* folder correctly:

- 1. Right-click the *Wheelsbridge* folder located in the *Program* (*x86*) folder and select *Properties*.
- 2. Select the **Security** tab.
- 3. In the upper combo text box select *User* and check if permission to read and write to the folder is granted (lower combo text box).
- 4. If not select *Edit*, select *User* in the upper combo text box.
- 5. Check permission to read and write in the lower combo text box.

- 6. Click OK.
- 7. The User can now read and write to the *Wheelsbridge* folder.