




Wheels Bridge



Nevus Analyzer TiVi102 User Manual

User Manual 1.2
Version 1.2
February 2019

PIONEERS IN TISSUE VIABILITY IMAGING

TiVi102 Nevus Analyzer

Dear Valued Customer!

Welcome to the WheelsBridge TiVi102 Nevus Analyzer toolbox intended for automatic analysis of naevi, moles and skin spots of different colours. The TiVi102 Nevus Analyzer was designed for observer-independent analysis of skin spots in association with evaluation of skin care products and for medical research applications.

The WheelsBridge TiVi 102 Nevus Analyzer utilizes a highly sensitive digital camera equipped with polarization filters making it possible to suppress surface reflections from the skin. The versatile system software – based on the MATLAB® high performance language for technical computing – allows for rapid and easy capturing and analysis of images. Among the many useful features of the TiVi 102 Nevus Analyzer software the following are of particular interest:

- Automatic capturing of photos.*
- Automatic tracking of skin objects by colour or intensity discrimination.*
- Removal of large and small objects from the nevus lesions map.*
- Individual display of objects (nevus lesions) by a scanning through procedure.*
- Individual determination of object area, border, colour and diameter.*
- Comparison of two naevi maps lesion by lesion.*
- Small objects removal from plot possible.*
- Automatic identification of spots with growing areas.*
- All data can be exported to ASCII-format spread sheets.*

We are convinced that the TiVi102 Nevus Analyzer will be a productive tool in the assessment and follow-up of human skin treatment with skin care products and in medical research.

Thank you for choosing the WheelsBridge TiVi102 nevus Analyzer.

WheelsBridge AB

TiVi102 Nevus Analyzer

Information in this document is subject to change without notice and does not represent a commitment on the part of WheelsBridge AB. No part of this manual may be reproduced or transmitted in any form or by any means, electronically or mechanically, including photocopying and recording, for any purpose without the written permission of WheelsBridge AB.

Any use of application of WheelsBridge products must be performed by properly trained and qualified personnel. WheelsBridge AB takes no responsibility for action taken on the basis of information derived from the use of the product described in this manual or from consequential damages in connection with furnishing, performance or use of this manual.

Licensee's limited rights are governed by the license agreement between WheelsBridge AB and MathWorks. Licensee may not modify or remove any license agreement file (MathWorks or third party) that is included with the MCR Libraries ("MCR Library License"). Licensee of this Application accepts the terms of the MCR Library License.

The *Nevus Analyzer* is not registered as a Medical Device. It is intended for research applications only.

WheelsBridge is a registered trade mark.

1. INTRODUCTION

Nevus (moles, birthmarks, brown spots) is the medical term for well-defined and benign lesions of the skin. Skin cancer – especially malignant melanoma – however may arise from pre-existing nevus lesions. It is therefore of interest to closely follow the development of nevus spots with respect to changes in colour, area, border and diameter over time and pay special attention to protecting sensitive skin from UV radiation and other skin damaging factors.

The *TiVi102 Nevus Analyzer* makes it possible to automatically track important nevus parameters such as area, border, colour and diameter and quantify how these parameters vary over time and over the skin surface. It is further possible to compare the change of individual naevi areas in two different photos and automatically detect spots with enlarged areas.

The intended use of the *TiVi102 Nevus Analyzer* is to analyse skin spots of different shapes and colour in experimental and research applications. It is not yet approved for the diagnosis and treatment of disease.

2. OPERATING PRINCIPLE

The basic operating principle of the *TiVi102 Nevus Analyzer* is probably best explained by a flow chart (Fig 1). After having captured a high-resolution photo (typically 18 Mpixels) of the central part of the skin of the back, this photo is uploaded to the *TiVi102 Nevus Analyzer*.




Fig 1. Basic operating principle Flow Chart.

1. Upload the photo to the *TiVi102 Nevus Analyzer*.
2. Zoom in and move the zoomed-in window over the photo until a representative nevus becomes visible.
3. Click on the nevus in the photo to display its colours against gray-scale background (the colour space around the centre colour can be selected by the user). Alternatively, the intensity of the nevus can be used for colour space discrimination.
4. Zoom out until the target size of the sub-photo is attained. All nevus lesions within the selected colour space become visible.

5. The individual nevus lesions are converted to individual objects. Their location can be displayed as a black and white map.
6. Too small and too large objects can be eliminated (the user sets the size elimination threshold).
7. The remaining individual objects (nevus lesions) can be displayed one-by-one. Their area, border, colour and diameter are calculated and displayed.
8. All calculated data can be exported to an Excel® spreadsheet.

3. PHOTO CAPTURING

When capturing photos for the *TiVi102 Nevus Analyzer*, use the highest camera resolution and place the camera about 1 m away from the object. This will give a field of view of about 430 mm x 290 mm and a resolution better than 0.1 mm per pixel. Set the filters to cross-polarized mode. Set the **AutoFocus** switch to AF and the **Stabilizer** to 'ON'. Use the zoom in feature to attain a field of view that captures a photo of the central part of the back (Fig. 2), in order to avoid shadows from the periphery and to get the highest resolution allowing also small skin objects to be included and analyzed.




Fig 2. Place the camera one meter away from the object and use the zoom in feature to attain a field of view that includes the central part of the skin of the back.

Open the *TiVi701Camera* window and set the **Size Reduction** to 'Minor' (under the **Size Reduction** pull-down menu), the **Select Photo Size** to 'Large Fine' and the **No of Photos** to '1' in accordance with Fig 3.




Fig. 3. Setting Camera parameters in the *TiVi701Camera* window.


Power the **Illuminator** and switch off the ambient light in the room. Capture a photo and store it in the selected folder under the selected name in the computer.

3. GETTING STARTED


The basic features of the *TiVi102 Nevus Analyzer* are probably best explained by way of an example. In the following example it is assumed that the photo *nevus-0001.jpg* has been captured by a TiVi camera system and stored in the *TiVi102 demonstration* folder. This photo displays the central part of a back with many nevus lesions of different sizes and colours.

1. Open the *Nevus Analyzer* toolbox from the **Toolbox** pull-down menu in the *TiVi700 Analyzer* window.

TiVi102 Nevus Analyzer



2. Click the **Load Photo** button and navigate to the *TiVi102 demonstration* folder and select the *nevus-0001.jpg*.




3. Move the **Zoom In** slider upwards to about 90% of its maximum value. The part of the photo displayed in the **Photo** panel is outlined in the **Image** panel photo.

TiVi102 Nevus Analyzer




4. Place the mouse pointer on the photo displayed in the **Photo** panel. Press the left mouse button and drag the mouse over the photo. Release the mouse button. Repeat the procedure until a representative nevus lesion appears in the zoomed-in photo.




TiVi102 Nevus Analyzer

- Click on the nevus lesion in the zoomed-in photo without moving the mouse. The selected nevus lesion will now be displayed against the surrounding skin displayed in gray colour scale.




- Move the **Zoom In** slider downwards to zoom out the photo to a size suitable for further analysis. The zoomed-out window can be moved over the original photo by placing the mouse pointer on the photo, pressing the left mouse button and dragging the mouse pointer to a new position. The site of the zoomed-out window is outlined in the original photo displayed in the **Image** panel.


TiVi102 Nevus Analyzer



7. Click the **Mask** radio-button to display the zoomed-out part of the photo in black and white scale.




- Click the **Mixed** radio-button and then the **Display Objects** button to open the **Objects** window.




- This sub-photo includes 261 objects. To eliminate the smallest objects, move the **Remove Small** slider to its rightmost position. The smallest objects are now deleted and only 9 objects remain in the zoomed-in window. The average area, border, colour and diameter of these 9 objects are displayed in the **Average** panel.

TiVi102 Nevus Analyzer




10. Click the **Mixed** radio-button to display the selected nevus lesion objects against a gray scale coloured skin background.




11. Click the **Enumerate** button to display the object number next to the object.

TiVi102 Nevus Analyzer



12. Click the **Mixed** radio-button in the **Display** panel. Then click the **Scan Through** button to display the different objects one-by-one. The parameters for the object displayed are shown in the **Actual Object** panel.



TiVi102 Nevus Analyzer

- Click the **Scan Through** button repeatedly to display the different objects and the associated parameters.
- Click the **Export Data** and select a name and folder to save the Excel file. Open the Excel file to display the data.


File Name NC-0002.jpg											
Red	Green	Blue	Expand	Mode	Rem Smal	Rem Large	Col1	Col2	Raw1	Raw2	
0,45898	0,225882	0,179922	0,1	Col. Mode	501	5000	1476	3429	511	1813	
Number	Avg Area	Avg Borde	Avg Colou	Avg Diameter							
9	1492,66	255,38	255,38	41,95							
Number	Area	Border	Diameter								
1	665	194,55	29,09								
2	1057	287,56	36,68								
3	2073	315,4	51,37								
4	1839	240,1	48,38								
5	2873	339,88	60,48								
6	2413	285,7	55,42								
7	1302	269,07	40,71								
8	685	210,24	29,53								
9	527	155,92	25,9								

- This completes the **GETTING STARTED** section.

4. COMPARISON OF NAEVI MAPS


In the following example it is explained how the areas of different naevi in two different photos can be compared and how an enlarging spot can be automatically identified.

1. Open the *Nevus Analyzer* toolbox from the **Toolbox** pull-down menu in the *TiVi700 Analyzer* window.
2. Click the **Load Photo** button, navigate to the *TiVi102demonstration* toolbox and double-click on the Act-0001 photo.




3. Click the **Load Ref Photo** and double-click the *Ref-0001* photo.

TiVi102 Nevus Analyzer




4. Place the mouse pointer on the darkest nevus click the left mouse button.




5. The **Mixed photo** is now displayed.

TiVi102 Nevus Analyzer




6. Click the **Mask** radio-button to display the black and white photo.




7. Drag the **Expand Color** slider somewhat to the right to expand the colour space and include more spots.

TiVi102 Nevus Analyzer




8. Click the **Show Actual Photo** radio-button to display the black and white actual photo.




9. Select **Synch Actual to Ref** in the **Select Standard Colour** pull-down menu.

TiVi102 Nevus Analyzer




10. Click the **Compare Images** button.




11. The reference spot (largest spot) is now indicated by a red square in the **Reference Photo**. Click on the corresponding spot in the **Actual** photo to synchronize the two photos.

TiVi102 Nevus Analyzer




12. Print '100' in the **Remove Small Objects** edit box to remove small objects.




13. Click the **Show Plot** button to display the *Plot 102* window. The **Diagram** in the *Plot102* window displays the relationship between spot sizes in the **Reference Photo** and the **Actual Photo**. Observe that one spot has a substantially higher value in the **Actual Photo** than in the **Reference Photo** as indicated by the red dot in the **Diagram**.


TiVi102 Nevus Analyzer



- Click on the red dot with the highest **Actual Area** value in the **Diagram**. The corresponding spot is now identified in the photo by circumscribing it with a red rectangle.



15. Check the **Mask** radio-button to display the corresponding black and white photo.



16. Check the **Refresh** button to hide the red rectangle circumscribing the actual spot. Check the **Actual Photo** radio-button and then the **Reference Photo** radio-button to verify that the actual spot is indeed larger in the **Actual** photo than in the **Reference** photo.
17. Click the **Export Data** button to export the results to an *Excel* file. Open the *Excel* file to display the data set.
18. This completes the **COMPARISON OF NAEVI MAPS** section.

5. DETAILED DESCRIPTION

The TiVi102 Nevus Analyzer main window




Photo Panel

1. **Photo text** – a text string that guides the user to the next step in the image processing process.
2. **Zoom In slider with edit box** – to zoom in on a specific part of the photo. The section displayed is outlined in the **Image Panel** photo.
3. **Clicking in the photo** – placing the mouse pointer somewhere in the photo, pressing the left mouse button and moving the mouse with the button pressed causes the zoomed in region to move over the original photo. Clicking on a spot inside the photo without moving the mouse generates a new photo with the colours representing the actual point displayed against a gray-scale background photo. The colour space – the colours displayed – is composed of colours around the colour of the point clicked on. The size of this colour space is determined by the setting of the **Expand Color** slider.
4. **Photo size** – displays the size of the zoomed-in photo.

5. **Date and Time** – displays the date and time of photo capturing.

Select Photos Panel

1. **Load Photo button** – to load a new photo from file.
2. **Load Ref Photo button** – to load a new reference photo from file.
3. **Show Actual Photo rb-button** – to display the actual photo.
4. **Show Ref Photo rb-button** – to display the reference photo.
5. **Display Panel: Photo, Mixed & Mask** – check the **Photo** radio-button to display the entire photo in colour. Check the **Mixed** radio-button to display the coloured nevus lesions against gray-scale background skin. Check the **Mask** radio-button to display the photos as black and white images.
6. **Actual Photo edit-box** – displays the name of the actual photo.
7. **Reference Photo edit-box** – displays the name of the reference photo.
8. **Display Objects button** – opens the **Object** window and transfers image information to this window.
9. **Clear All button** – deletes photos and resets the system to start mode.
10. **Display Colour Space button** – opens the **Colour Space** window.
11. **Reset button** – resets the settings but keeps the photos uploaded and displayed.
12. **Coordinates Panel** – displays the coordinates of the zoomed-in photo.

Image Panel


1. **Image text** – displays the selected mode: **Colour Mode** indicates that the **Colour Space** is selected with the centre colour widened by the settings of the **Expand Colour** slider. In this mode a dark nevus lesion can contain a core area outside the selected **Colour Space**. **Intensity Mode** (default) indicates that the **Intensity** is selected with the centre intensity increased by the setting of the **Expand Color** slider. In this mode all pixels of an intensity lower than the centre intensity increased by the **Expand Color** slider setting are included in the nevus lesions identified, i.e. a dark nevus lesion will include its centre darker colours in the mask.
2. **Centre edit-box** – displays the centre colour of the **Colour Space**.
3. **High edit-box** – displays the brightest colour included in the **Colour Space**.

4. **Low edit-box** – displays the darkest colour included in the **Colour Space**.
5. **Size** – displays the size of the photo in pixels.

Select Colour Panel

1. **Red slider & edit-box** – sets the red centre colour (automatically updated when a pixel is pointed on in the photo in the **Photo Panel**).
2. **Green slider & edit-box** – sets the green centre colour (automatically updated when a pixel is pointed on in the photo in the **Photo Panel**).
3. **Blue slider & edit-box** – sets the blue centre colour (automatically updated when a pixel is pointed on in the photo in the **Photo Panel**).
4. **Expand Colour slider & edit-box** - sets the extension of the **Colour Space** around the centre colour (**Colour Mode**) or intensity threshold level above the centre intensity (**Intensity Mode**).
5. **Select standard Colour pull-down menu** – sets predetermined centre colours and width of **Colour Space** (and updates the sliders and edit boxes accordingly). *By Mouse Pointer* – the **Color Space** is set by pointing at a point in the photo and clicking the left mouse button. *Light Nevus* – The RGB-values of the centre point are set to [0.5 0.3 0.2] and the **Expand Color** slider value to 0.1. *Dark Nevus*– The RGB-values of the centre point are set to [0.25 0.15 0.1] and the **Expand Color** slider value to 0.1. *Synchronize* – the **Colors Space** of the photo in focus is synchronized (made equal to that of) with the hidden photo. After having clicked on a nevus to define a colour space, two more items are added to the pull-down-menu if both the actual and the reference photo are uploaded. **Synch Act to Ref** synchronizes the colour space of the actual photo to that of the reference photo, while **Synch Ref to Act** synchronizes the colour space of the reference photo to that of the actual photo. After synchronizing the colour spaces of the two photos the **Compare Images** button is enabled.
6. **Compare Images button** – to open the *Compare 102* window.
7. **Int Mode button** – this toggle button sets the system to operate in **Colour Mode** or **Intensity Mode**.
8. **Load Settings** – loads saved parameter settings from file.
9. **Save Settings** – saves parameters to file.
10. **Close** – closes the *TiVi102 Nevus Analyzer Toolbox*.

The Colour Space window




Limits Panel

1. **Red Low Limit** – displays the low red limit of the **Colour Space**.
2. **Red High Limit** – displays the high red limit of the **Colour Space**.
3. **Green Low Limit** – displays the low green limit of the **Colour Space**.
4. **Green High Limit** – displays the high green limit of the **Colour Space**.
5. **Blue Low Limit** – displays the low blue limit of the **Colour Space**.
6. **Blue High Limit** - displays the high blue limit of the **Colour Space**.

Select Colour Space Panel

1. **ROI rb-button** – displays the results in the **RGB Colour Space**.
2. **HSV rb-button** – displays the results in the **HSV Colour Space**.
3. **Close-button** – closes the **Colour Space** window.
4. **Diagram rotation** – to rotate the diagram by pointing on it with the mouse pointer and pressing the left mouse button while moving the mouse.

The Object window



Object Panel

1. **Photo text** – a text string that guides the user to the next step in the image processing process.
2. **Photo** – displays the zoomed-in-photo imported from the *TiVi102 Nevus Analyzer* main window.


Data Panel

1. **File Name** – displays the filename.
2. **Actual Object Panel** – displays number, area, border, colour and diameter of the object (nevus lesion) selected.
3. **Average Panel** – displays number of nevus lesions and average area, border, colour and diameter of the nevus lesions in the zoomed in photo.

Control Panel


1. **Remove Small Object slider & edit-box** – removes small objects from the zoomed-in photo.
2. **Remove Large Object slider & edit-box** – removes large objects from the zoomed-in photo.
3. **No.Bin edit-box** – to set the number of bins in the histogram.
4. **Histogram Type pull-down menu** – to select the type of histogram to be displayed (area, border or diameter).
5. **Show Histogram button** – toggles between **Show histogram / Show Image** to display histogram or image.
6. **Display Object button** – toggles between **Display Object / Display Map** to display individual objects (single nevus lesions) or image map.
7. **Number of Objects text** – displays the number of objects in the zoomed-in photo.
8. **Size of Map text** – size of map displayed.
9. **Mode text** – indicates mode of operation – **Intensity Mode** or **Colour Mode**.
10. **Display Panel: Photo, Mixed & Mask** – displays the zoomed-in photo as pure colour photo, coloured nevus lesions against gray background skin or black and white object mask.
11. **Indicator Settings Panel** – to set the size of numbers and their distance to the objects in the black and white image.
12. **Enumerate button** – toggles between **Enumerate / No Numbers** to display or hide object numbers in the black and white image.
13. **Export Data button** – to export data to Excel® file.
14. **Set Object Size slider and edit-box** – to set the size of the individual object photo. With the slider set to zero the boundaries of the largest object fills the entire photo.
15. **Scan Through button and edit box** – scans through and displays the individual objects while updating the **Actual Object** panel parameters accordingly.
16. **Close button** – to close the **Object** window.

TiVi102 Nevus Analyzer



Object window with a single object (nevus lesion) displayed.

The Compare 102 window



Actual Photo Panel

1. Displays **File Name** and **Number of Objects** in **Actual Photo** .

Reference Photo Panel

1. Displays **File Name** and **Number of Objects** in **Reference Photo** .

Controls Panel

1. **Remove Small Objects slider** - to remove small objects from the photos .
2. **Remove Small Objects edit** - to remove small objects from the photos .
3. **Enumerate** – to display number of the objects in the reference photo.
4. **Show Plot** – to open the *Plot102* window.
5. **Close** – to close the **Compare 102** window.

Display Panel

1. **Photo radio-button** – to display the full colour photos.
2. **Mixed radio-button** – to display the spots in colour space colours superimposed on the full photo displayed in gray-scale colour.
3. **Mask radio-button** – to display the spots in white colour on black background.

The Plot 102 window




Diagram Panel

1. Displays areas of spots in the actual photo versus the corresponding areas of the spots in the reference photo. Actual photo spot sizes are displayed as red dots if the corresponding actual spot – reference spot areas ratio exceeds the warning limit as set by the **Warning slider**. Spots along the y-axis (spots existing only in the actual photo) are indicated by yellow dots while spots along the x-axis (spots existing only in the reference photo) are indicated by green dots. The remaining spots are indicated by blue dots.

2. Any dot in the **Diagram** clicked on by the mouse pointer automatically identified the corresponding spot in the reference photo by circumventing it by a red square in the photo.

Photo Panel

1. Displays the photo, mixed photo or the spots (white on black background). A red square is automatically drawn around the spot corresponding to the dot clicked on by the mouse in the **Diagram**.

Display Panel

1. Displays filename of the actual and reference photo displayed as well as the number of objects in the corresponding photos.
2. **Photo radio-button** – to display the full colour photos.
3. **Mixed radio-button** – to display the spots of selected colour space colours against a gray-scale version of the photo.
4. **Mask radio-button** – to display the spots in the actual or reference photo in white colour against a black background.
5. **Actual photo** – to display the **Actual Photo**.
6. **Reference photo** – to display the **Reference Photo**.
7. **Export Data** – to export data to an *Excel* spread-sheet.
8. **Refresh** – to remove the red square circumventing the selected spot from the photo.
9. **Close** – to close the *Plot 102* window.